PoliLed PoliScaler

FARO

L'ODONTOTECNOLOGIA

PoliLed

LAMPADA POLIMERIZZANTE A LED

Covered by Design Patent

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GB

Poliled is a curing light for dental units that Faro has realised with LED technology.

Its spare design is attractive and ergonomic.

The control electronics, achieved with SMT technology, are fully integrated in the handpiece.

With limited dimensions (l = 165mm, max Ø 27mm), **Poliled** is made of aluminium alloy and weighs only 100 g.

Poliled is equipped with a single button that is used to select, activate and halt the emission when desired.

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GB

An acoustic signal is emitted at the start, mid-point and end of the selected cycle.

Two different curing cycles can be performed with **Poliled**.

By pressing the button normally the power cycle is selected, which enables the emission of maximum power for the entire duration (20").

Holding the button down for 2" selects the soft cycle, which makes the emission gradually ramp up to the maximum (in 10"), then maintains it for another 10" (total duration 20").

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GB

Pressing the button during the curing interrupts the emission.

Poliled uses a high-efficiency 5W LED that is able to supply a power of over 1500 mW/cm², measured at the outlet of the optical fibre.

The optical fibre has a silicone rubber tip that prevents even accidental contact with the surface to be cured.

The optical fibre and the tip can be sterilised in an autoclave at 134°C.

Poliled comes with a dedicated power supply cord to which it is fastened by a connector that allows it to be

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GB TECHNICAL SPECIFICATIONS	
Application:	Polymerisation of light-hardening dental materials
Power supply:	From 17 to 26 Vac 50Hz From 19 to 33 Vdc
Consumption:	11.5 VA
Light source:	One 5W LED with optical concentrator
Wavelength:	440 – 480 nm
Typical intensity:	1500 mW/cm ² 700 mW/cm ² (in the Compact version)
Optical fibre:	Ø 8 autoclavable at 134°C (500 cycles)
Instrument size:	Length: 165 mm Max. diam.: 27 mm
Cord:	Detachable by means of connector. Integral with handpiece (in the Compact version)
Weight:	100 g
Material:	Aluminium
Control electronics:	Fully integrated in the instrument
Service life:	1.8 million 20-second cycles
Applied standards:	EN 60601-1 EN 60601-1-2
Working times:	<u>Standard cycle:</u> 20 sec. at maximum power <u>Ramp cycle:</u> 10 sec. from zero to maximum power + 10 sec. at maximum power <u>Acoustic signals:</u> At start and every 10 sec

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GB

detached when needed (for maintenance or other operations).

A more economical version exists, with a power of 700 mW/cm² and the cord integral with the polymeriser, called **Poliled compact**.

Poliled is also made with a special attachment that makes it possible to connect it to the unified cord which is able to power both the **Poliled** and the **PoliScaler** ablator.

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GB

PoliScaler is a new ultrasound scaler manufactured by FARO. In modern dentistry, ultrasound scaler are now part of the instrumentarium that every professional should have as a technological aid for carrying out their activity as efficiently as possible. Its ergonomic shape, reduced size and light weight (max. diam. 17 mm – length without tip 104 mm – weight 41 g) makes it possible to use **PoliScaler** while retaining perfect visibility of the operating field. The smooth surface of the handpiece facilitates cleaning prior to sterilisation in an autoclave at 134°C. For sterilisation a special case is provided that in addition to the handpiece can hold the tips used along with their supports. **PoliScaler** comes with three different types of tips (D2-D3-D4) which cover most operating requirements. Upon request another two tip models can be supplied (D1-D5) for more specific applications. Replacing the tips is simplified by the special support with which each tip is equipped. Using the supplied torque wrench ensures optimum

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GB

fastening of the tip on the handpiece. **Poliscaler** is available in the version for the dental unit and the portable version (**P-Mobile**). In the dental unit version, the electronic control, mounted on an aluminium support, can be easily fixed in the unit thanks to its reduced dimensions (h = 34mm, l = 60mm, d = 54mm). The ultrasound intensity can be adjusted to optimize its effectiveness by using the special potentiometer in the dental unit version, and by using the appropriate buttons in the **P-Mobile** version. Operation of the **PoliScaler** and activation of the cooling spray are controlled by the pedal of the dental unit (in the built-in version) and by a small dedicated pedal in the portable version. To protect the patient's health, **Poliscaler** is provided with an anti-suction valve that prevents the contaminated water from pushing up into the cord by blocking it in the handpiece, which can be sterilised, thereby avoiding the risk of contamination of the next patient. Another important feature is the possibility of mounting on the cord of the

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GB

PoliScaler the **Poliled** curing light in the special version, in the Compact and normal models (articles 650800 and 650810). This solution allows both instruments to be always available and immediately operative by attaching them, without distinction, to the cord and selecting the button of the instrument desired.

Specific functions of the tips

Type D1:

- For the zone between the tooth surface and the gum tissue.
- Supragingival, subgingival and interdental removal of tartar with slow vertical movement.
- The Type D1 tip must be used with an angle of less than 90° with respect to the tooth.
- The Type D1 tip is normally used with low power.

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GB

Type D2:

- For scaling in general.
- Supragingival, subgingival and interdental removal of tartar and stains.
- The Type D2 tip is normally used with medium power.
- For correct use of the Type D2 tip, hold it at a 30° tangent to the tooth.
- The Type D2 tip must not be used with an angle of 90° with respect to the tooth.

Type D3:

- For treatment of periodontal pockets and interdental spaces, especially of teeth with more than one root.
- Normally Type D3 is used with low power.
- Insert the Type D3 tip into the periodontal pockets to measure their depth (max. 8mm).

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GB TECHNICAL SPECIFICATIONS	
Device classification:	Class 2
Applied classification rule:	Typo B
Membership class 93/42 ECC:	II a
Operation mode:	Continuous
Delivery pressure of water:	From 1 to 5 bar
Water consumption:	100÷120 cc./min.
Supply voltage:	24 Vac, 50/60 Hz
Absorption:	24 VA
Frequency of ultrasound unit:	24-32 kHz
Weight of handpiece:	41 g
Water temperature :	Mains temperature

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GB
<ul style="list-style-type: none"> - Use the Type D3 tip, with slow tangential movements, along the whole tooth for 30-60 sec. - Do not press on the surface of the tooth. <p>Type D4:</p> <ul style="list-style-type: none"> - For the zone between the tooth surface and the gum tissue. - The Type D4 tip is longer than the Type D1 tip. - The Type D4 tip must not be used at 90° with respect to the tooth surface. - Usually Type D4 tip is used with low power.

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GB TECHNICAL SPECIFICATIONS	
Device classification:	Class 1
Applied classification rule:	Typo B
Membership class 93/42 ECC:	II a
Operation mode:	Continuous
Delivery pressure of water:	From 1 to 5 bar
Water consumption:	100÷120 cc./min.
Power supply voltage:	Input AC 100-240V, 50/60Hz 1,0A Output DC 24V 1,7A
Absorption:	24 VA
Frequency of the ultrasound unit:	24-32 kHz
Weight of handpiece:	41 g
*Water temperature :	Mains temperature
Degree of protection against water penetration:	Control IP 40 Transformer IP 20

* water loading tube from mains complete with connectors included

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GB

Type D5:

- For scaling of very hard tartar and removal of cement.
- The Type D5 tip can be used vertically on the surface of the tooth.
- Usually the Type D5 tip is used with medium power.
- For the type of applications of the Type D5 tip, it is particularly important to check for proper delivery of water in order to ensure cooling of the tip.

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FARO S.p.A. via Faro, 15 20060 Ornago (Milano) ITALIA - Tel. +39 039 68781 - Fax + 39 039 6010540